

WOLF-PAK™



*HuberPRO*®

# Safety Huber Needle Winged Infusion Set



- **SIMPLE**
- **SAFE**
- **SECURE**

## FEATURES

---

- Easy to use
- Low profile
- Audible and visible confirmation of safety feature
- Latex free
- Non-DEHP
- Available with or without the needle free access Y-Site

## BENEFITS

---

- Minimizes risk of needlestick injury
- Low profile allows for maximum patient comfort
- Color coded wing for ease of needle gauge identification



WOLF MEDICAL SUPPLY, INC.

EXCLUSIVE PRODUCTS AND OUTSTANDING SERVICE  
TO ORDER: 800.335.9653 FAX: 954.835.0301  
[www.wolfmed.com](http://www.wolfmed.com)



# Safety Huber Needle Winged Infusion Set

**Needle safety doesn't have to be complicated.**



**1**

Hold as shown and remove Sheath.



**2**

Insert per protocol.



**3**

Flip guard forward to lay flat against the patient. Secure per protocol.



**4**

Push guard forward with index finger to capture the needle as it exits from the patient. An audible "click" will notify you that the needle is captured.

## DESCRIPTION

## ORDER NUMBER

## SHELF BOX QTY/CASE QTY

Micro-bore 19G x 1" w/o y-site	HP1910	25/300
Micro-bore 19G x 1" w/y-site	HP1910Y	25/300
Micro-bore 19G x 3/4" w/o y-site	HP1975	25/300
Micro-bore 19G x 3/4" w/y-site	HP1975Y	25/300
Micro-bore 20G x 1" w/o y-site	HP2010	25/300
Micro-bore 20G x 1" w/y-site	HP2010Y	25/300
Micro-bore 20G x 3/4" w/o y-site	HP2075	25/300
Micro-bore 20G x 3/4" w/y-site	HP2075Y	25/300
Micro-bore 22G x 1" w/o y-site	HP2210	25/300
Micro-bore 22G x 1" w/y-site	HP2210Y	25/300
Micro-bore 22G x 3/4" w/o y-site	HP2275	25/300
Micro-bore 22G x 3/4" w/y-site	HP2275Y	25/300
Standard-bore 19G x 1" w/y-site	HPS1910Y	25/300
Standard-bore 19G x 3/4" w/y-site	HPS1975Y	25/300
Standard-bore 20G x 1" w/y-site	HPS2010Y	25/300
Standard-bore 20G x 3/4" w/y-site	HPS2075Y	25/300



WOLF MEDICAL SUPPLY, INC.

**EXCLUSIVE PRODUCTS AND OUTSTANDING SERVICE**  
**TO ORDER: 800.335.9653 FAX: 954.835.0301**  
[www.wolfmed.com](http://www.wolfmed.com)