



WOLF-PAK® BIOSITE PLUS® NEEDLELESS VALVES

NEUTRAL SYSTEM

Our Wolf-Pak® Biosite Plus® needleless valves feature neutral displacement connector technology, which reduces the risk of bacterial contamination and promotes patient safety.

WHY WOLF-PAK® BIOSITE PLUS®?

- A single protocol is to be used with all patients, making it a seamless transition to purchase and implement clinically
- For use on all vascular catheters for blood draws or administration of IV medications
- A mechanically and microbiologically closed system reduces the risk of infection
- MRI compatible
- Improves hospital efficiencies and helps lower IV therapy costs

Promotes aseptic technique. Contains no latex and is non-DEHP.
Pressure rated to psi levels listed below.



OPAQUE

The opaque valve is pressure rated for use on PICC lines.

In-Stock Available Sizes (all packaged 100/Box):

BI8000
BI8500

Biosite Plus® Clear Needleless Valve
Biosite Plus® Opaque Needleless Valve

(60 psi)
(300 psi)



WOLF MEDICAL SUPPLY, INC.

EXCLUSIVE PRODUCTS AND OUTSTANDING SERVICE
TO ORDER: 800.335.9653 | FAX: 954.835.0301

www.wolfmed.com



WOLF-PAK® BIOSITE PLUS® NEEDLELESS VALVES

NEUTRAL SYSTEM

ANTIMICROBIAL BARRIER

Elastomeric silicone design minimizes the migration of microorganisms through it.

MINIMAL DEAD SPACE

Minimal priming volume of 0.04 mL accommodates lower flush volumes.

POLYCARBONATE SHAFT

Sterile, transparent spike permits continuous and laminar flow, preventing occlusion.

FLUID PATHWAY

Zero fluid displacement eliminates connection and disconnection reflux.

SEPTUM

Flat, smooth, and swabbable external surface allows for safe clockwise and anticlockwise cleaning.

EXTERNAL BODY

Smooth and ergonomic polycarbonate protection layer provides visual confirmation of flush after use with medications or blood.

Easy grip design minimizes slips and touch contamination.



CLEAR

The clear housing allows for visualizing the internal fluid path when flushing the connector.



WOLF MEDICAL SUPPLY, INC.

EXCLUSIVE PRODUCTS AND OUTSTANDING SERVICE
TO ORDER: 800.335.9653 | FAX: 954.835.0301
www.wolfmed.com