



An easy way to select the perfect flow rate for patients:

FREEDOM60[®] uses flow rate tubing to control the rate of drug delivery. These economical 20 piece boxes contain a variety of our most popular flow rate tubing sets.

There's no need to stock excessive inventory of all different FREEDOM60[®] flow rate tubing sets to be sure you'll have exactly what new patients need. Starter Kits make it convenient to begin therapy using one flow rate, and then increase or decrease drug flow rate by switching to a different tubing set as appropriate for the patient.

L20KT Low Flow Rate Starter Kit

Contains:
(3) F30
(2) F45
(8) F60
(4) F120
(3) F180

H20KT High Flow Rate Starter Kit

Contains:
(2) F275
(5) F600
(5) F900
(4) F1200
(4) F2400



FREEDOM60[®] Syringe Infusion System

- ▶ Advanced Technology – No batteries or electricity required!
- ▶ Low Cost per Use – Uses standard 60cc syringes
- ▶ Reliable – No maintenance or calibration cycles required
- ▶ Flow rates available for a variety of time frames, ranging from a few minutes to 5 days.
- ▶ Automatically reduces flow rates as sites saturate
- ▶ No runaway, freeflow, bolus or overdose
- ▶ Constant Safe Pressure – Holds pressure during/after infusion



F60_FL_tubing_starter_kits_03-1003vB

FREEDOM60[®] is CE Approved.

RMS Medical Products fully complies with the ISO 13485 Quality Management System.