



WOLF MEDICAL SUPPLY, INC.

Contact us today for more information
800.335.9653 or orders@wolfmed.com
www.wolfmed.com

WOLF-PAK® Smartsite®

Needle-Free Systems

Wolf Medical Supply offers a full line of Needle-Free I.V. disposables to support the needs of your facility and patients.



- ◆ Metal-free
- ◆ MRI compatible
- ◆ DEHP-free



WOLF-PAK®



Safety is our first priority...

- ◆ No needles, no caps
- ◆ Airtight seal between the male luer and the top of the blue piston of the valve
- ◆ Smooth top is easy to swab
- ◆ Reduces the risk of needlestick injuries
- ◆ Compatible with standard luer connection for delivering safe and convenient IV therapy with complete confidence
- ◆ Straight through fluid path helps improve the ability to flush
- ◆ Highly chemical compatible
- ◆ Cost-effective and lipid resistant
- ◆ Can be used with low pressure power injectors up to 325 psi and maximum flow rate of 10mL/second