

**IN STOCK  
ORDER TODAY!  
800.335.9653**



**WOLF MEDICAL SUPPLY, INC.**

**ALL SIZES  
AVAILABLE  
ORDERS@WOLFMED.COM**

**WOLF-PAK®**



**TPN BAGS**

**Features:**

- Reinforced hanger
- Detachable transfer sets
- Pre-attached, non-reopening clamp
- Made from 100% EVA material
- Sterile fluid pathway
- Easy Filling
- Compounding & Gravity filling


**Compounder bags compatible with:**

- Hospira Nutrimix
- Baxter Automix
- Healthtech Ad Comp
- Secure Medical Auto Comp



**TPNCP1L 1000ML COMPOUNDER  
TPNCP2L 2000ML COMPOUNDER  
TPNCP3L 3000ML COMPOUNDER**

**50/CASE**



**TPNGR1L 1000ML GRAVITY 3 LEG  
TPNGR2L 2000ML GRAVITY 3 LEG  
TPNGR3L 3000ML GRAVITY 3 LEG**

**50/CASE**



**WOLF MEDICAL SUPPLY, INC.**

**www.wolfmed.com**