




800.335.9653  
 orders@wolfmed.com  
 www.wolfmed.com


# A Size for Every Situation



Grip-Lok® Universal Securement Device	Actual Length	Minimum Tubing Size	Maximum Tubing Size
---------------------------------------	---------------	---------------------	---------------------


Featuring 3 sizes to secure virtually all lines, tubes and catheters.





	
1/16" 1.5mm	3/16" 4.5mm




	
1/8" 3mm	5/16" 8mm




	
1/4" 6mm	1/2" 13mm



Grip-Lok H® Hydrocolloid Securement Device	Actual Length	Minimum Tubing Size	Maximum Tubing Size
--	---------------	---------------------	---------------------

Featuring a hydrocolloid adhesive that is **gentle** on patients with fragile skin.



	
1/32" 1mm	1/8" 3mm



	
1/16" 1.5mm	3/16" 4.5mm



## WOLF MEDICAL SUPPLY, INC.